

AF6 Medium, modified
(Watanabe *et al.* 2000)

This is a good general freshwater medium that supports the growth of a wide diversity of algae. The original medium was designed for growing *Colacium* (Euglenophyceae) (Kato 1982). It has been modified at the NIES Culture Collection of Microalgae and Protozoa (Watanabe *et al.* 2000) by adding MES buffer, removing CaCO₃ and substituting a different trace metal solution.

Prepare the stock solutions. To 950 mL of dH₂O, first dissolve the MES and citrate, then add the stock solutions and bring the final volume up to 1 liter. The pH is adjusted to 6.6. Autoclave.

Component	Stock Solution	Quantity	Molar Concentration in Final Medium
MES buffer		400 mg	2.05 x 10 ⁻³ M
NaNO ₃	140 g L ⁻¹ dH ₂ O	1 mL	1.65 x 10 ⁻³ M
NH ₄ NO ₃	22 g L ⁻¹ dH ₂ O	1 mL	2.75 x 10 ⁻⁴ M
MgSO ₄ 7H ₂ O	30 g L ⁻¹ dH ₂ O	1 mL	1.22 x 10 ⁻⁴ M
K ₂ HPO ₄	5 g L ⁻¹ dH ₂ O	1 mL	2.87 x 10 ⁻⁵ M
KH ₂ PO ₄	10 g L ⁻¹ dH ₂ O	1 mL	7.35 x 10 ⁻⁵ M
CaCl ₂ 2H ₂ O	10 g L ⁻¹ dH ₂ O	1 mL	6.80 x 10 ⁻⁵ M
Fe-citrate	2 g L ⁻¹ dH ₂ O	1 mL	8.17 x 10 ⁻⁶ M
Citric acid	2 g L ⁻¹ dH ₂ O	1 mL	1.04 x 10 ⁻⁵ M
trace metals solution	(see below)	1 mL	---
Vitamin solution	(see below)	1 mL	---

Trace Metals Solution

This trace metals solution is a modified version of PIV trace metals solution (Provasoli and Pintner 1960). Prepare the primary stock solutions. Into 950 mL of dH₂O, pH EDTA to 8.0-8.4 to assure it fully dissolves, add Fe and fully, and finally add 1 mL of each primary stock.

Component	Primary Stock Solution	Quantity	Molar Concentration in Final Medium
Na ₂ EDTA 2H ₂ O	---	5.000 g	1.34 x 10 ⁻⁵ M
FeCl ₃ 6H ₂ O	---	0.98 g	3.63 x 10 ⁻⁶ M
MnCl ₂ 4H ₂ O	180 g L ⁻¹ dH ₂ O	1 mL	9.10 x 10 ⁻⁷ M
ZnSO ₄ 7H ₂ O	110 g L ⁻¹ dH ₂ O	1 mL	3.83 x 10 ⁻⁷ M
CoCl ₂ 6H ₂ O	20.0 g L ⁻¹ dH ₂ O	1 mL	8.41 x 10 ⁻⁸ M
Na ₂ MoO ₄ 2H ₂ O	12.5 g L ⁻¹ dH ₂ O	1 mL	5.17 x 10 ⁻⁸ M

AF6 Vitamin Solution

First, prepare the three primary stocks. Into 950 mL of dH₂O, add and dissolve the thiamine, and then add 1 mL from each of the two primary stocks. Filter sterilize. Refrigerate or freeze.

Component	Primary Stock Solution	Quantity	Molar Concentration in Final Medium
thiamine (vit. B ₁)	---	10 mg	2.96 x 10 ⁻⁸ M
biotin (vit. H)	2.0 g L ⁻¹ dH ₂ O	1 mL	8.19 x 10 ⁻⁹ M
cyanocobalamin (vit. B ₁₂)	1.0 g L ⁻¹ dH ₂ O	1 mL	7.38 x 10 ⁻¹⁰ M
pyridoxine (vit. B ₆)	1.0 g L ⁻¹ dH ₂ O	1 mL	5.91 x 10 ⁻⁹ M

Kato, S. 1982. Laboratory culture and morphology of *Colacium vesiculosum* Ehrb. (Euglenophyceae). *Jap. J. Phycol.* **30**: 63-67. (in Japanese with English summary).

Provasoli, L. and Pintner, I.J. 1960. Artificial media for fresh-water algae: problems and suggestions. pp. 84-96. In Tyron, C.A. Jr. and Hartman, R.T. (eds.) *The Ecology of Algae*. Special Publ. 2, Pymatuning Laboratory of Field Biology, Univ. Pittsburgh, PA.

Watanabe, M.M., Kawachi, M., Hiroki, M. and Kasai, F. (Eds.) 2000. NIES Collection List of Strains. 6th Ed. NIES, Japan. 159 pp.

As of 2017 the NCMA has used standard f/2 vitamins without deleterious effects.

f/2 Vitamin Solution

(Guillard and Ryther 1962, Guillard 1975)

First, prepare primary stock solutions. To prepare final vitamin solution, begin with 950 mL of dH₂O, dissolve the thiamine, add the amounts of the primary stocks as indicated in the quantity column below, and bring final volume to 1 liter with dH₂O. At the NCMA we autoclave to sterilize. Store in refrigerator or freezer.

Component	Primary Stock Solution	Quantity	Molar Concentration in Final Medium
thiamine · HCl (vit. B ₁)	---	200 mg	2.96 x 10 ⁻⁷ M
biotin (vit. H)	0.1 g/ L dH ₂ O	10 mL	2.05 x 10 ⁻⁹ M
cyanocobalamin (vit. B ₁₂)	1.0 g/ L dH ₂ O	1 mL	3.69 x 10 ⁻¹⁰ M

Guillard, R.R.L. 1975. Culture of phytoplankton for feeding marine invertebrates. pp 26-60. In Smith W.L. and Chanley M.H (Eds.) *Culture of Marine Invertebrate Animals*. Plenum Press, New York, USA.

Guillard, R.R.L. and Ryther, J.H. 1962. Studies of marine planktonic diatoms. I. *Cyclotella nana* Hustedt and *Detonula confervacea* Cleve. *Can. J. Microbiol.* **8**: 229-239.